

qLabs[®] ElectroMeter Plus

for Professional Point of Care



Lightweight, portable, compact and easy to use, the qLabs platform features:

- Single-use and cost effective disposable strips
- Lightweight and compact handheld meter
- High precision with lab quality: 5% CV
- Fingertick blood: less than 10 μ l
- Rapid testing: 2-7 minutes
- Wireless communications for easy data uploading

Features

- Reliable laboratory results in less than 2 minutes for PT and 7 minutes for APTT
- 10 µL of fingertip whole blood, fresh venous whole blood or citrate whole blood
- On board quality control
- External liquid QC (LQC) available
- Built-in temperature controlled testing
- Built-in barcode scanner
- No calibration or instrument manipulation required



Current Test Menu

- PT
- aPTT
- PT/aPTT Combo

Interface

- Display 2.1 inch 128×96 LCD
- Barcode scanner Laser class 2 according to EN 60825-1:2007

Software

- Language Multi-language menu

Memory Capacity

- Patient records Stores up to 200 patient records and 12 QC results

Performance

- Measurable Range INR between 0.5 to 8.0 APTT from 20 sec to 420 sec

Specifications

- Dimensions 135mm x 65mm x 34.5mm (L× W × H)
- Weight 145 g (Without Batteries)
- Operating conditions Temperature: 10 °C - 35 °C
- Relative humidity: 10 % - 90 %
- Power supply 100V-240V AC
- 4 X AA batteries



Vet Model Also Available



Plot No 31 B kauvery Street, Ram Nagar
South, Madipakkam, Chennai - 600091
TamilNadu, india.
74490 92777, 74490 93777, 74490 95777

Email : info@quicklab.in,
Web : www.quicklab.in

Our Authorized Dealer

